Industrial Weighing Systems

9 Richmond St. Picton, ON Canada K0K 2T0

Ph: 613-786-0016 Cell: 613-921-0397 Fax: 613-476-5293

E-mail info@iwsystems.ca Website: www.iwsystems.ca

Industrial Weighing Systems assumes no responsibility for the accuracy of the calibration instructions extracted from manufactureres documentation.

Contact us for complete manuals we may have on file info@iwsystems.ca

**Industrial weighing systems provides repairs to all insturments and load cells.
Calibration services available.**

Should you require additional information please do not hesitate to contact IWSystems
Check out our products at www.iswystems.ca

ELECTRONIC SCALE

BILANCIA ELETTRONICA

Series
Serie

NL

OPERATING MANUAL
MANUALE OPERATIVO

LANGUAGES

LINGUE

ENGLISHpage 3

ITALIANOpag.22

INDEX

1. Precautionary measures	page 4
2. Installation	page 4
2.1 Installation of the scale	page 4
2.2 Installation of the windscreen	page 5
2.3 Preparation.....	page 5
3. Technical Specifications	page 6
3.1 Power supply.....	page 6
3.2 List of divisions	page 6
3.3 Dimensions	page 7
4. Display description	page 7
5. Keyboard functions	page 8
6. Setup of the scale	page 9
7. Functioning	page 16
8. User calibration	page 17
9. RS-232 interface	page 18
9.1 Transmission string	page 18
9.2 Examples of data transmission.....	page 19
10. Weight conversion and total capacity	page 20
10.1 Total capacities for each unit of measure	page 20
11. Error messages	page 21
DECLARATION OF CONFORMITY	page 41
WARRANTY and AUTHORISED ASSISTANCE CENTRE	page 41


1. Precautionary measures

- Do not wash under water; if any water comes in contact with the scale dry off immediately with a dry cloth.
- Install the scale in a dry place and away from heat sources.
- Put the load on the plate without causing knocks. Do not exceed the maximum indicated capacity.
- Never use solvents or acetone to clean the scale.
- Remove the batteries if the scale is not used for long periods of time.

2. Installation


2.1 Installation of the scale

- Remove the scale from the box handling it carefully.
- Put the scale on a stable surface. The weighing area must remain clean and dry. Install the scale away from extreme heat or cold. Protect the scale from vibrations, dust and air currents.
- Put the plate on its appropriate support without causing knocks, making sure not to press the plate.
- Remove the transportation protection which is on the right side of the scale. Keep this part for eventual future transportations.
- Adjust the feet of the scale in order to improve its stability (check the level).
- Connect the scale to the power supply mains, using the fitted power adapter.
- Press the switch-on key on the right side of the scale.
- Wait for the scale to finish executing the Self Test and for the zero to appear on the display.

2.2 Installation of the windscreen

The windscreen is made up of 5 separate pieces. The 4 sliding sides together form a square, and the cover can be put over it as an extra protection.

2.3 Preparation

- Make sure that the level is in the centre.
- It's better to use a separate plug in order to avoid interferences with other electrical equipment.
- Do not put weights on the plate in the moment the scale is turned on.
- Put the weight in the centre of the plate. The weight should not exceed the maximum capacity of the scale.
- To obtain a better accuracy one should turn on the scale 15-20 minutes before using it.
- If the "+ -" indication appears on the upper right-hand corner of the display; it is advisable to change immediately the batteries.
- Operating temperature: 5 ~ 40°C.

3. Technical Specifications

3.1 Power supply

The scale is supplied with a 9V/300mA AC/DC Power adapter.

Low battery warning:

When the “+ -” symbol appears on the right corner of the display, one should change the battery.

Note: it's possible that when the low battery symbol appears on the display, the scale does not weight accurately and the zero is unstable.

3.2 List of divisions


Unit of measure	NL600H	NL150A	NL300A	NL1500A	NL3000A
g	600×0.01	150×0.005	300×0.01	15000×0.05	3000×0.1
ct	3000×0.05	750×0.05	1500×0.05	7500×0.5	15000×0.5
lb	1.3225×0.00005	0.3304×0.00002	0.66×0.00005	3.304×0.0002	6.6×0.0005
oz	21.15×0.0005	5.29×0.0002	10.575×0.005	52.9×0.002	105.75×0.005
dr	338.6×0.01	84.6×0.005	169.2×0.01	846×0.05	1692×0.1
gn	9258×0.2	2313×0.1	4626×0.2	23130×1	46260×2
ozt	19.29×0.0005	4.82×0.0002	9.64×0.005	48.2×0.002	96.4×0.005
dwt	385.8×0.01	96.4×0.005	192.8×0.01	964×0.05	1928×0.1
MM	160×0.05	40×0.002	80×0.005	400×0.02	800×0.05
tl.J	16×0.0005	4×0.0002	8×0.0005	40×0.002	80×0.005
tl.T	16×0.0005	4×0.0002	8×0.0005	40×0.002	80×0.005
tl.H	16×0.0005	4×0.0002	8×0.0005	40×0.002	80×0.005
t	51.4×0.001	12.85×0.0005	25.7×0.001	128.5×0.005	257×0.01

3.3 Dimensions:

Body	172(W)×230(D)×60(H)mm
Plate	Thickness : 7mm Diameter : 114mm
Windscreen	150(W)×162(D)×82(H)mm
Net weight	1030 g

4. Display description

Display size: 84mm x 24mm

5. Keyboard Functions

Cali. Scale calibration key.

Unit Unit of measure selection key.

Zero Display zeroing key


Tare Tare key.

Mode Serves to select the functioning mode (weighing, counting, percentile).

Set Confirmation key.

↑ Selection key.

Enter Confirmation key.


6. Setup of the scale

Through the Setup environment, it's possible to configure the scale according to the user needs. It's possible to select the specific unit of measure (up to 13 different units), check the speed of the display update, the RS-232 interface and the configuration of the zero tracking.

To enter the setup of the scale:

Keep the **Cali.** key pressed and simultaneously press the **Enter** key. The display will show the first parameter, "SPEEd".

1. Press **Enter** to enter in the desired parameter; press it again to save the eventual modifications made and pass to the following parameter.
2. Press the **Set** key to choose the desired function.
3. Press the **Tare** key to exit the Setup by saving any changes made.
4. Press the **Zero** key to exit the Setup, without saving any changes made.


“Speed” parameter: Display speed update

Press the ↑ key to choose the available speeds of the display update. There are 5 different choices: “**SPd 1**” lower updating speed, “**SPd 5**” greater updating speed.

Press **Enter** to confirm the choice and pass to the following parameter, or press **Set** to not save and pass to the following parameter.

(The default value is “**SPd 2**”).


“Zero” parameter: Zero tracking range

Press the \uparrow key to choose the available ranges of the zero tracking. There are 5 different choices: “ZEr 1” lower range, “ZEr 5” upper range (see the table at the end of the page).

Press **Enter** to confirm the choice and pass to the following parameter, or press **Set** to not save and pass to the following parameter.

(The default value is “ZEr 1”).

Zero tracking values:

	NL150A, NL300A, NL1500A, NL3000A series	NL600H series
“Zer 1”	0.8d	1.6d
“Zer 2”	1.5d	2.5d
“Zer 3”	2.5d	3.5d
“Zer 4”	3.5d	4.5d
“Zer 5”	4.5d	5.5d


“b-Lgt” parameter: Not used in this application

“Unit” parameter: Selectable unit of measure in weighing

Press the ↑ key to choose the available units of measure. In order to be able to select the unit of measure in weighing, select “ON”, and to disable the unit of measure select “OFF”.

Press **Enter** to confirm the choice and pass to the following unit of measure, or press **Set** to not save and pass to the following parameter.

(The default value is “on”).


“InitU” parameter: Unit of measure selection upon start-up.

Press **Unit** to scroll the available units of measure.

Press **Enter** to confirm the choice and pass to the following parameter, or press **Set** to not save and pass to the following parameter.

(The default value is “g”.)

“trn” parameter: RS232 transmission mode

Press the ↑ key to choose the transmission mode:

“Stb”= transmission upon weight stability.


“Etb” = transmission of the weight upon pressing **Enter** (with stable weight).

“SEr”= continuous weight transmission.

“CLOSE” = disabled transmission.

Press **Enter** to confirm the choice and pass to the following parameter, or press **Set** to not save and pass to the following parameter.

(The default value is “CLOSE”).


“rAtE” parameter: RS232 Baud Rate

Press the ↑ key to set the desired Baud Rate: “1200” “2400” “4800” “9600”.

Press **Enter** to confirm the choice and pass to the following parameter, or press **Set** to not save and pass to the following parameter.

(The default value is **2400**).

7. Functioning

Turn on the scale by pressing the on key on the left side. The scale will be ready for weighing, by using the unit of measure used previously. Press the **Mode** key to choose the weighing, counting or percentile (%) function.

1 Weighing

1.1 SELECTION OF THE UNIT OF MEASURE.

Press the **Unit** key to choose the unit of measure: the display will show the selected value. Up to 13 units of measure can be enabled.

1.2 TARE FUNCTION

Put the container on the scale and, when the weight is stable, press **Tare**. The display will show the zero weight (net) and will store the tare weight. To cancel the tare function press **Tare** without any weight on the plate.

1.3 ZERO FUNCTION

If the scale with an unloaded plate shows the gross Zero, press **Zero** to reset it. The maximum clearable field is $\pm 4\%$ of the capacity.

2 Counting

The scale will count the pieces weighing a preset number of samples on the plate. If necessary, put a container on the plate and press **Tare** before starting to weigh the samples.

2.1 REFERENCE OPERATION

- a) Press **Mode** once: the display shows "SA 10pcs".
- b) Select the number of pieces with the \uparrow key (10, 20, 50, 100, 200, 500, 1000) and press **Enter** to confirm (the display shows "PCS").
- c). Put the pieces to be counted on the scale, wait for the weight stability indication and press **Enter**.
- d) If the sampling has been made correctly, the display will show the number of pieces selected in point b.

Note: If the weight used for the sampling is less than 80% of the minimum division of the scale, the display will show: "**-CSL- PCS**". Press **Enter** to return to weighing or press **Set** to carry out a new sampling.

3 Percentile mode

3.1 SAMPLING WITH THE WEIGHT

- a) Press **Mode** twice: the display shows “-PRE-%”.
- b) Put the weight on the scale, wait for the stability indication and press **Enter** to confirm.
- c) If the sampling has been executed correctly the display will show “100%”; now the percentile of the weights will be relative to that used for the sampling.

Note: If the weight used for the sampling is less than 0.1% of the total capacity; the display will show: “-PSL-%”. Press **Enter** to return to weighing or press **Set** to execute a new sampling.

3.2 SAMPLING WITHOUT WEIGHTS

Press the **Set** key. The display will show “-PRE-%”, then press the **Unit** key.

Manually enter the value of the sample weight: with the \uparrow key the blinking digit increases and with the **Enter** key one confirms and passes to the following digit; after the confirmation of the last digit, the percentile function will enable with the sampling just made.

Note:

Once the sampling procedure in the counting or percentile mode (%) has been made, press the **Mode** key to switch the display between WEIGHT >> PIECES >> PERCENTILE. The sampling data will remain stored.

8. User calibration

- 1) Press the **Cali.** key until the scale shows “-000-”
- 2) Press **Enter** to acquire the Zero of the scale; after a few instants the display will show the sample weight for the calibration. By pressing the \uparrow key it's possible to select one of the proposed sample weights: 1/6, 1/3, 1/2, 2/3, 5/6, 1 of the total capacity.
- 3) Put the selected calibration weight on the plate, wait for the weight stability and confirm with **Enter**. After a few instants the scale will automatically go in the weighing mode.

9. RS232 interface

The scale has a standard RS232 interface for the data output.

Connector: 9 pin
Pin 2: TRANSMISSION
Pin 5: GND

The output format may not be modified. Each transmitted data is codified in 10 bits:

bit 1, bit 2, bit 3, bit 4, bit 5, bit 6, bit 7, bit 8, bit 9, bit 10.

bit 1: start bit
bit 2~bit 9: data
bit 10: stop bit

9.1 Transmission string

The transmission string is made up of 18 bytes:

Meaning	Length in bytes	Contents
Functioning mode	2	"WT" Weighing. "CT" Counting "PC" Percentile
Weight status	2	"OL" weight in over load. "ST" stable weight. "US" unstable weight.
Sign	1	"+" "-"
Weight	7	The decimal point is included.
Percentile unit of measure.	4	Spaces in counting mode
String end	2	CR (0D hexadecimal), LF (0A hexadecimal)

The transmission is made up of a start bit, 8 data bits, one stop bit.

9.2 Examples of data transmission:

1. Net stable weight of 38.25 g:

W	T	S	T	+				3	8	.	2	5				g	0D	0A
---	---	---	---	---	--	--	--	---	---	---	---	---	--	--	--	---	----	----

2. Unstable net weight of 300 ct:

W	T	U	S	+					3	0	0				c	t	0D	0A
---	---	---	---	---	--	--	--	--	---	---	---	--	--	--	---	---	----	----

3. Stable net weight of - 60.0 tl.H:

W	T	S	T	-				6	0	.	0	t	l	.	H		0D	0A
---	---	---	---	---	--	--	--	---	---	---	---	---	---	---	---	--	----	----

4. +60 pieces, with stable net weight:

C	T	S	T	+						6	0						0D	0A
---	---	---	---	---	--	--	--	--	--	---	---	--	--	--	--	--	----	----

5. + 25,4%, with stable net value:

P	C	S	T	+				2	5	.	4						0D	0A
---	---	---	---	---	--	--	--	---	---	---	---	--	--	--	--	--	----	----

10. Weight conversions

1 ct (Carats) = 0.1999694 g
 1 lb (Pounds) = 453.59237 g
 1 oz (ounce) = 28.349523 g
 1 dr (DRAM) = 1.7718451 g
 1 gn (GRAIN) (UK) = 0.0647989 g
 1 ozt (TROY OUNCE) = 31.103476 g
 1 dwt (PENNYWEIGHT) = 1.555174 g
 1 MM (MOMME, Japan) = 3.749996 g
 1 tl. J (JEWELRY TAEI, Hong Kong) = 37.42900 g
 1 tl. T (TAEI, Taiwan) = 37.49995 g
 1 tl. H (TAEI, Singapore) = 37.799375 g
 1 t (TOLA, India) = 11.663804 g

10.1 Total capacities for each unit of measure.

	NL600H	NL150A	NL300A	NL1500A	NL3000A
g	600.09	150.045	300.09	1500.45	3000.9
ct	3000.45	750.45	1500.45	7504.5	15004.5
lb	1.32295	0.33058	0.66045	3.3058	6.6045
oz	21.545	5.2918	10.5795	52.918	105.795
dr	338.69	84.645	169.29	846.45	1692.9
gn	9259.8	2313.9	4627.8	23139	46278
ozt	19.2945	4.8218	9.6445	48.218	96.445
dwt	385.89	96.445	192.89	964.45	1928.9
MM	160.045	40.018	80.045	400.18	800.45
T1.J	16.045	4.0018	8.0045	40.018	80.045
t1.T	16.045	4.0018	8.0045	40.018	80.045
t1.H	16.045	4.0018	8.0045	40.018	80.045
t	51.409	12.0045	25.709	120.045	257.09

11. Error messages

- ✧ “**0-Err**” : Means that when the scale is turned on there is a weight greater than 20% of the capacity; remove the weight and turn the scale back on.
- ✧ “**-OL-**” : Means that a weight greater than the value of the capacity + 9 divisions has been put on the scale. The scale will emit an acoustic signal.
- ✧ “**-LO-**” : Insufficient power supply voltage; the indication will appear with a weight equal to 0.
- ✧ “**UNSTA**” : After the scale is turned on, an unstable weight has been detected.
- ✧ “**-Adc-**” : The AD converter is in overflow condition; an acoustic alarm is emitted.
- ✧ “**-PSL-**” : In the percentile mode, an insufficient sample weight has been used.
- ✧ “**-SLAC-**” : Means that, in the counting mode, the sample weight is less than 50% of the scale division.
- ✧ “**-CSL-**” : means that, in the counting mode, the sample weight is less than 80% of the scale division.

RECYCLING INSTRUCTIONS

The crossed-out wheeled bin on the product means that at the product end of life, it must be taken to separate collection or to the reseller when a new equivalent type of equipment is purchased. The adequate differentiated refuse collection in having the product recycled, helps to avoid possible negative effects on the environment and health and supports the recycling of the materials of which the equipment is made. The unlawful disposal of the product by the user will entail fines foreseen by the current regulations.

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This device conforms to the essential standards and norms relative to the applicable European regulations. The Declaration of conformity is available in the web site www.scalehouse.it.

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WARRANTY GARANZIA

The TWO-YEAR warranty period begins on the day the instrument is delivered. It includes spare parts and labour for repairs at no charge if the INSTRUMENTS ARE RETURNED prepaid to the DEALER'S PLACE OF BUSINESS. Warranty covers all defects NOT attributable to the Customer (such as improper use) and NOT caused during transport.

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If on site service is requested (or necessary), for any reason, where the instrument is used, the Customer will pay for all of the service technician's costs: travel time and expenses plus room and board (if any).

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The customer pays for shipping costs (both ways), if the instrument is shipped to the DEALER or manufacturer for repair.

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The WARRANTY is VOIDED if faults occur due to work done by unauthorised personnel or due to connections to equipment installed by others or incorrect connection to the power supply.

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This warranty DOES NOT provide for any compensation for losses or damages, direct or indirect, incurred by the Customer due to complete or partial failure of instruments or systems sold, even during the warranty period.

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